

Med Pack Contents with Supply Locations

Compartment Designations:



It is suggested that you lay all supplies on a table and group them by assigned compartments

Then, methodically put those supplies in each location.

Remember to squeeze air out of all plastic bags.

Med Pack Contents with Supply Locations

Compartment		Qty
A	Small Pocket (far left) see Diagram	
	Gloves -SM, Exam	4 pr.
	Gloves - MED, Exam	4 pr.
B	Zipper Pocket (in front of C) see Diagram	
	LIST - <i>Med Pack Contents with Supply Locations</i>	1
C	Medium Pocket (left, split) see Diagram	
Front	Eye & Wound Wash (4 oz /15ml) Exp. Date	2
	Hand Sanitizer (Alcohol Foam - Purell eq) 2oz	1
	Kelly Forceps 5 1/2" (Locking Hemostat)	1
	Splinter Forceps, Fine point, Serrated 3 1/2" ("Tweezers")	1
Back	Eye & Wound Wash (4 oz /15ml) Exp. Date	1
	Hand Sanitizer (Alcohol Foam - Purell eq) 2oz	1
	Thermometer, Veridian 2-Second Digital (w/ Instructions)	1 Batteries
	Penlight w/ Pupil Gauge, reuseable (2 AAA batteries)	1 Batteries
D	Zipper Pocket (in front of E) see Diagram	
	Combat Application Tourniquet (CAT)	1
E	Large Pocket (center) see Diagram	
	Olaes Vacuum Sealed Bandages - 6" (Dressing, Field) (hemorrhage control bands for moderate bleed)	2
F	Rollout Pouch (4 Compartments, center)	
1	ABD Trauma Dressing, 10x30 or 12x30	1
	ABD Trauma Pads, 5x9 or 8x10	2
	Mask, N95 (flat fold)	3
	SAM Splints, 4.25x36 (foam)	1
	Bio-Hazard Bags, red	2
	Eyewear, Protective - required over Rx Glasses	1
2	Triangular Bandages	2
	Shears (EMS Trama, 7 1/4")	1

F Rollout Pouch (4 Compartments, center)			
3	Blanket, Mylar Survival	2	
4	ABD Trauma Pads, 5x9 or 8x10	4	
Wet	Alcohol Wipes	20	
Wet	Povo-Iodine Prep Pads	10	
Wet	BZK Antiseptic Towelette	20	
Wet	Triple Antibiotic Ointment Packs	10	
Dry	Eye Pad, Oval Dressing	6	
Dry	Gauze Pads, 4x4 or 3x3 (Can be Non-Sterile)	20	Stack
Dry	Telfa Pads, 2x3" or 3x4" (Non Adherent)	5	
Dry	Gauze, stretch 3", 4 yds.	3	rolls
Dry	Bandage, 1"x3" Fabric	10	
Dry	Bandage, X-Large Fabric, J&J Band-Aid	10	
Dry	Butterfly Bandage (Wound Closure)	10	
Dry	Steri-Strip (Wound Closure)	6	
Dry	Gauze, Petrolatum or Occlusive, 3x9	1	

G Zipper Pocket (in front of H) see Diagram			
***** E M P T Y *****			

H Medium Pocket (right, split) see Diagram			
Front	Coban 3", 5 yds. (self-adhesive)	2	rolls
	1" Medical Paper Tape	1	roll
Back	Coban 3", 5 yds. (self-adhesive)	2	rolls
	2" Medical Paper Tape	1	roll

I Small Pocket (far right) see Diagram			
	Gloves - LRG, Exam	4	pr.
	Gloves - XLG, Exam	4	pr.